



Panel	Preanesthetic	Renal	Diagnosis Plus	Equine	Large Animal	Liver Plus	Avian & Reptile	Diabetes	Diagnosis-II	Critical Care
ALB	X	X	X	X	X	X	X	X	X	X
ALP	X		X	X	X	X			X	X
ALT	X		X			X			X	X
AMY			X						X	
AST				X	X	X	X			
BA						X	X			
BUN	X	X	X	X	X	X			X	X
CA		X	X	X	X		X		X	X
CHOL						X			X	
CL		X			X		X			X
CPK				X	X		X			X
CREA	X	X	X	X					X	X
FRU								X		
GGT				X	X	X			X	
GLU	X		X	X	X		X	X	X	X
K		X	X	X	X		X			X
LAC										X
LIPA									X	
MG					X					
NA		X	X	X	X		X			X
PHOS		X	X		X		X		X	
TBIL			X	X		X			X	
tCO2		X		X						X
TP	X		X	X	X	X	X	X	X	X
TRIG										
UA							X			
A/G#	X		X	X	X		X	X	X	X
AGap#		X								X
B/C#	X	X	X	X					X	X
C-Ca#		X	X	X	X		X		X	X
GLOB#	X		X	X	X	X	X	X	X	X
Na/K#		X	X	X	X		X			X
UREA#	X	X	X	X	X	X			X	X

Einzelparameter
BUN
CREA
LAC
NH ₃
TBIL
UPC

Dualparameter
ALB + TP
ALP + ALT
AST + CPK
BA + CHOL
CA + PHOS
GGT + CPK
LIPA + AMY
TRIG + CHOL

Immuno-parameter
cCRP
SAA
TT4

#berechnete Parameter