

# scil Vet abc Plus<sup>+</sup>

## Referenzwerte



Parameter	konv. Einheit	Hund	Katze	Pferd	Kaninchen	Schwein	Schaf	Rind	Maus	Ratte	Frettchen
<b>WBC</b>	10 <sup>3</sup> /mm <sup>3</sup>	6.0 - 12.0	5 - 11	5 - 10	6.3 - 10.0	3.5 - 10.0	3.5 - 10.0	7 - 10	3 - 15	3 - 15	2.5 - 5.5
<b>RBC</b>	10 <sup>6</sup> /mm <sup>3</sup>	5.5 - 8.5	5 - 10	5.5 - 11.0	5.2 - 6.8	3.8 - 5.8	3.8 - 5.8	5.5 - 8.0	5 - 12	5 - 12	6.5 - 12.0
<b>HGB</b>	g/dl	15 - 20	8 - 17	9 - 18	11.5 - 15.1	11.0 - 16.5	11.0 - 16.5	10 - 13	11.1 - 18.0	11.1 - 18.0	15.0 - 18.0
<b>HCT</b>	%	44 - 57	27 - 47	30 - 50	36 - 47	35 - 50	35 - 50	30 - 40	36 - 52	36 - 52	40.0 - 60.0
<b>MCV</b>	µm <sup>3</sup>	60 - 77	40 - 55	37 - 55	65 - 76	80 - 97	80 - 97	45 - 60	44 - 69	44 - 69	45.0 - 65.0
<b>MCH</b>	pg	17 - 26	13 - 17	13 - 19	21.1 - 24.5	26.5 - 33.5	26.5 - 33.5	16 - 22	12.0 - 24.5	12.0 - 24.5	15.0 - 20.0
<b>MCHC</b>	g/dl	31 - 38	31 - 36	31 - 39	29.5 - 33.9	31.5 - 35	31.5 - 35.0	28 - 36	21.6 - 42.0	21.6 - 42.0	30.0 - 34.0
<b>RDW</b>	%	14 - 17	17 - 22	14 - 20							14 - 17
<b>PLT</b>	10 <sup>3</sup> /mm <sup>3</sup>	200 - 460	180 - 430	100 - 400	250 - 610	150 - 390	150 - 390	175 - 600	140 - 600	140 - 600	300 - 900
<b>MPV</b>	µm <sup>3</sup>	6.7 - 11.1	6.5 - 15	6.7 - 11.1							6.7 - 11.1
<b>LYM #</b>	10 <sup>3</sup> /mm <sup>3</sup>	1.0 - 3.6	1.0 - 4.0	1.5 - 4.0							0.3 - 1.3
<b>MON #</b>	10 <sup>3</sup> /mm <sup>3</sup>	0.0 - 0.5	0.0 - 0.5	0.0 - 0.4							0.0 - 0.2
<b>GRA #</b>	10 <sup>3</sup> /mm <sup>3</sup>	3.0 - 10.0	3.0 - 12.0	3.0 - 8.0							0.4 - 2.0
<b>EOS #</b>	10 <sup>3</sup> /mm <sup>3</sup>	0.0 - 0.6	0.0 - 0.6	0.0 - 0.4							0.0 - 0.27



**scil animal care company GmbH**  
**Dina-Weissmann-Allee 6**  
**68519 Viernheim**  
**Tel.: +49 (0) 6204 78 90 221**  
**Fax: +49 (0) 6204 78 90 200**  
**Email: info-def@scilvet.com**  
**www.scilvet.com**

# scil Vet abc Plus<sup>+</sup>

## Referenzwerte



Parameter	SI Einheit	Hund	Katze	Pferd	Kaninchen	Schwein	Schaf	Rind	Maus	Ratte	Frettchen
<b>WBC</b>	10 <sup>9</sup> /l	6 - 12	5 - 11	5 - 10	6.3 - 10.0	3.5 - 10.0	3.5 - 10.0	7 - 10	3 - 15	3 - 15	2.5 - 5.5
<b>RBC</b>	10 <sup>12</sup> /l	5.5 - 8.5	5 - 10	5.5 - 11.0	5.2 - 6.8	3.8 - 5.8	3.8 - 5.8	5.5 - 8.0	5 - 12	5 - 12	6.5 - 12.0
<b>HGB</b>	g/l	150 - 200	80 - 170	90 - 180	115 - 151	110 - 165	110 - 165	100 - 130	111 - 180	111 - 180	150 - 180
<b>HCT</b>	l/l	0.44 - 0.57	0.27 - 0.47	0.30 - 0.50	0.36 - 0.47	0.35 - 0.50	0.35 - 0.50	0.30 - 0.40	0.36 - 0.52	0.36 - 0.52	0.4 - 0.6
<b>MCV</b>	fl	60 - 77	40 - 55	37 - 55	65 - 76	80 - 97	80 - 97	45 - 60	44 - 69	44 - 69	45.0 - 65.0
<b>MCH</b>	pg	17 - 26	13 - 17	13 - 19	21.1 - 24.5	26.5 - 33.5	26.5 - 33.5	16 - 22	12.0 - 24.5	12 - 24.5	15.0 - 20.0
<b>MCHC</b>	g/l	310 - 380	310 - 360	310 - 390	295 - 339	315 - 350	315 - 350	280 - 360	216 - 420	216 - 420	300 - 340
<b>RDW</b>	%	14 - 17	17 - 22	14 - 20							14 - 17
<b>PLT</b>	10 <sup>9</sup> /l	200 - 460	180 - 430	100 - 400	250 - 610	150 - 390	150 - 390	175 - 600	140 - 600	140 - 600	300 - 900
<b>MPV</b>	fl	6.7 - 11.1	6.5 - 15.0	6.7 - 11.1							6.7 - 11.1
<b>LYM #</b>	10 <sup>9</sup> /l	1.0 - 3.6	1.0 - 4.0	1.5 - 4.0							0.3 - 1.3
<b>MON #</b>	10 <sup>9</sup> /l	0.0 - 0.5	0.0 - 0.5	0.0 - 0.4							0.0 - 0.2
<b>GRA #</b>	10 <sup>9</sup> /l	3.0 - 10.0	3.0 - 12.0	3.0 - 8.0							0.4 - 2.0
<b>EOS #</b>	10 <sup>9</sup> /l	0.0 - 0.6	0.0 - 0.6	0.0 - 0.4							0.0 - 0.27



**scil animal care company GmbH**  
**Dina-Weissmann-Allee 6**  
**68519 Viernheim**  
**Tel.: +49 (0) 6204 78 90 221**  
**Fax: +49 (0) 6204 78 90 200**  
**Email: info-def@scilvet.com**  
**www.scilvet.com**

# scil Vet abc Plus<sup>+</sup>

## Referenzwerte



Parameter	SI Einheit [mmol/l]	Hund	Katze	Pferd	Kaninchen	Schwein	Schaf	Rind	Maus	Ratte	Frettchen
<b>WBC</b>	10 <sup>9</sup> /l	6 - 12	5 - 11	5 - 10	6.3 - 10.0	3.5 - 10.0	3.5 - 10.0	7 - 10	3 - 15	3 - 15	2.5 - 5.5
<b>RBC</b>	10 <sup>12</sup> /l	5.5 - 8.5	5.0 - 10.0	5.5 - 11.0	5.2 - 6.8	3.8 - 5.8	3.8 - 5.8	5.5 - 8.0	5 - 12	5 - 12	6.5 - 12.0
<b>HGB</b>	mmol/l	9.3 - 12.4	4.97 - 10.56	5.59 - 11.18	7.14 - 9.38	6.83 - 10.25	6.83 - 10.25	6.21 - 8.07	6.89 - 11.18	6.89 - 11.18	9.3 - 11.2
<b>HCT</b>	l/l	0.44 - 0.57	0.27 - 0.47	0.30 - 0.50	0.36 - 0.47	0.35 - 0.50	0.35 - 0.50	0.30 - 0.40	0.36 - 0.52	0.36 - 0.52	0.4 - 0.6
<b>MCV</b>	fl	60 - 77	40 - 55	37 - 55	65 - 76	80 - 97	80 - 97	45 - 60	44 - 69	44 - 69	45.0 - 65.0
<b>MCH</b>	fmol	1.06 - 1.43	0.81 - 1.06	0.81 - 1.18	1.31 - 1.52	1.65 - 2.08	1.65 - 2.08	0.99 - 1.37	0.75 - 1.52	0.75 - 1.52	0.93 - 1.24
<b>MCHC</b>	mmol/l	19.25 - 23.60	19.25 - 22.36	19.25 - 24.22	18.32 - 21.05	19.56 - 21.73	19.56 - 21.73	17.39 - 22.36	13.41 - 26.08	13.41 - 26.08	18.6 - 21.1
<b>RDW</b>	%	14 - 17	17 - 22	14 - 20							14 - 17
<b>PLT</b>	10 <sup>9</sup> /l	200 - 460	180 - 430	100 - 400	250 - 610	150 - 390	150 - 390	175 - 600	140 - 600	140 - 600	300 - 900
<b>MPV</b>	fl	6.7 - 11.1	6.5 - 15.0	6.7 - 11.1							6.7 - 11.1
<b>LYM #</b>	10 <sup>9</sup> /l	1.0 - 3.6	1.0 - 4.0	1.5 - 4.0							0.3 - 1.3
<b>MON #</b>	10 <sup>9</sup> /l	0.0 - 0.50	0.0 - 0.5	0.0 - 0.4							0.0 - 0.2
<b>GRA #</b>	10 <sup>9</sup> /l	3.0 - 10.0	3.0 - 12.0	3.0 - 8.0							0.4 - 2.0
<b>EOS #</b>	10 <sup>9</sup> /l	0.0 - 0.6	0.0 - 0.6	0.0 - 0.4							0.0 - 0.27



**scil animal care company GmbH**  
**Dina-Weissmann-Allee 6**  
**68519 Viernheim**  
**Tel.: +49 (0) 6204 78 90 221**  
**Fax: +49 (0) 6204 78 90 200**  
**Email: [info-def@scilvet.com](mailto:info-def@scilvet.com)**  
**[www.scilvet.com](http://www.scilvet.com)**